

BP-201M

Wrist-type Fully Automatic Blood Pressure Monitor




One Touch Operation

- **Reliable**
- **Clinically Accurate**
- **120 Memories**

BP-201M

Wrist-type Fully Automatic Blood Pressure Monitor



SYSTOLIC

"M" Button

DIASTOLIC

"Power" Button

PULSE RATE

- ACCURATE MEASUREMENT
- QUICK READ
- QUIET INFLATION METHOD
- 2 X AAA BATTERIES
- AUTOMATIC POWER OFF


ISO9001
ISO13485 CE 0197

Product Specification


Measurement Method	Oscillometric Method	
Measurement Range	Pressure	0mmHg ~ 300mmHg
	Pulse	30 to 180 Beats/Minute
Measurement Accuracy	Pressure	± 3mmHg
	Pulse	± 5%
Pressurization	Automatic Pressurization	
Memory	120 Memories	
Function	Low Battery Detection	
	Automatic Power-Off	
Power Source	2 x AAA Batteries	
Unit Weight	Approx. 117 g (4.13oz.) (Excluding Battery)	
Unit Dimensions	Approx. 76 x 67 x 28mm (3" x 2.64" x 1.12")(L x W x H)	

Easy to use


Step 1
Wrap cuff around wrist.



Step 2
Press "Power" button to start testing.



Step 3
Press "M" button to view previously recorded measurements.



Owner's manual included